



We can provide interpreters for a variety of languages, information in larger print or other formats (e.g. audio) - please call us on 01932 723553.

You can use Text Relay to turn telephone communications into text.  
Use the Relay UK app on your phone, tablet or PC.  
You can also use Relay UK via Minicom or Uniphone.

اگر نیاز به ترجمہ دارید، لطفاً با شماره 01932 723553 تماس بگیرید.

ਜੇ ਤੁਹਾਨੂੰ ਤਰਜਮੇ ਦੀ ਲੋੜ ਹੈ ਤਾਂ ਕਿਰਪਾ ਕਰਕੇ ਇਸ ਨੰਬਰ 'ਤੇ فون کریں: 01932 723553

اگر آپ اس کا اردو زبان میں ترجمہ چاہتے ہیں، تو براہ کرم اس فون نمبر 01932 723553 پر رابطہ کریں

Se precisa de uma tradução por favor contacte: 01932 723553

আপনার অনুবাদের দরকার হলে এখানে যোগাযোগ করুন : 01932 723553

यदि आपको अनुवाद की ज़रूरत है तो कृपया इस नंबर पर फोन करें: 01932 723553

Jeżeli chcemy, aby te informacje w innym języku, proszę zadzwonić 01932 723553

**Ashford Hospital**  
London Road  
Ashford, Middlesex  
TW15 3AA  
Tel: **01784 884488**

**St. Peter's Hospital**  
Guildford Road  
Chertsey, Surrey  
KT16 0PZ.  
Tel: **01932 872000**

Website: [www.ashfordstpeters.nhs.uk](http://www.ashfordstpeters.nhs.uk)

# Omron ECG Monitor

## Clinical Measurements Department



## What is an Omron?

An Omron is a portable electrocardiogram (ECG) monitor that is used to record your heart rate and rhythm. It is a handheld device which can be used to record your symptoms. These symptoms need to last at least a few seconds from onset to allow you time to hold the monitor to your bare skin.

## What happens when I come for my appointment?

***Please ask a member of staff if you require a chaperone during your appointment.***

Before being shown how to use the device, an ECG will be performed. This test requires you to remove your clothing to the waist and have electrodes (stickers) attached to your chest, wrists, and ankles. This procedure takes approximately five minutes and is painless. You will then be shown how to use the device. This will take about five minutes to explain, after which you may then get dressed.

## What happens when I have an episode?

You will need to carry the monitor with you at all times. Please do not allow the monitor to get wet. When you need to make a recording, you need to **switch on** the monitor and place your index finger closely on the two finger electrodes at one end of the monitor. You then place the chest electrode (on the other end of the monitor) on **bare skin** about 5 cm below the left nipple.

When you are ready to take a measurement, you press the **START** button on the front of the monitor while maintaining the same posture. The monitor will beep every second for about

30 seconds – you will need to keep as still as you can and not move until the measurement is complete in order to create a clear reading.

The recording ends automatically, which is indicated by 4 rapid beeps. You may then switch off the monitor. Please ignore any messages that the device may display.

## What happens after the test?

You will be provided with, and ask to sign, a waiver to return the monitor. If the monitor is lost or damaged, you may be charged for a replacement. You will need to return the device to **Clinical Measurements Department in St Peters Hospital** where you picked up the monitor on the date that you are given when you attend your appointment. Please call the department on **01932 722530** if you have any concerns.

The results will be sent to the Consultant who requested the test. The results will be reviewed, and a letter will be sent to your General Practitioner (GP) advising them of the results. You will receive a copy of this letter. This may take 2-4 weeks.

---

## Further Information

We endeavor to provide an excellent service at all times, but should you have any concerns please, in the first instance, raise these with the Manager on duty. If they cannot resolve your concern, please contact our Patient Advice and Liaison Service (PALS) on 01932 723553 or email [asp-tr.patient.advice@nhs.net](mailto:asp-tr.patient.advice@nhs.net). If you remain concerned, PALS can also advise upon how to make a formal complaint.

Author: C Camargo

Department: Clinical Measurements

Version: 3

Published: April 2024

Review: April 2026

---